

## Choosing the right size has an impact on the result of the treatment with FroggyMouth!

The lips should not be stretched by the froggy Mouth in the relaxed, closed state. So the distance between the corners of the mouth should be slightly larger than the selected model.

Technical explanation:


Choosing a smaller size allows the lips not to be stretched. Stretching would result in non-physiological proprioception. The goal of the Froggy Mouth is to inhibit the motor and sensory activity of the facial nerve that stimulates the lips.

Vertically, there are no differences:

The vertical gap of the Froggy Mouth does not differ depending on the size of the appliance. The gap corresponds to the physiological opening of the mouth and does not allow the patient to create lip closure, which is the first step of the suck-swallow sequence.

As a rule, the choice of the FroggyMouth size depends on the clinical experience of the physician/therapist, taking into account the information listed above.

However, if in doubt, we recommend that you choose a smaller device. A smaller device is never less efficient, it's all about comfort.



Should any uncertainties arise, we are happy to clarify them:

Melanie Haslinger  
info@logicana.at  
+43 676 7053463

